

## Product information: SPY555-golgi (SC208)

Live Cell Fluorogenic Probe for the Golgi apparatus


### Introduction

SPY555-Golgi is a bright, far red & non toxic live cell golgi probe based on our SPY™ dyes series. Its optimized structure allows quick labeling of the golgi apparatus in live cells with high specificity and low background. SPY555-Golgi stains the golgi apparatus in live cells without the need for genetic manipulation or overexpression of fluorescent proteins. Its absorbance and emission spectra are similar to Cy5. SPY555-Golgi enables multicolor imaging with SPY505, SPY595, SPY650, SPY700 or GFP. SPY555-Golgi can be imaged with standard Cy3/TMR settings. It can be used for, confocal, SIM or STED imaging in living cells and tissue. Contains 1 vial of SPY555-Golgi (lyophilized).



### Probe Properties

<b>Absorbance maximum <math>\lambda_{abs}</math></b>	555 nm
<b>Fluorescence maximum <math>\lambda_{fl}</math></b>	580nm
<b>Works on fixed cells?</b>	yes
<b>Probe quantity</b>	100 stainings*
<b>Fluorescence lifetime</b>	2.4 ns
<b>STED depletion wavelength</b>	660 or 775 nm
<b>Shipping</b>	room temperature
<b>Storage</b>	-20°C



### Storage & Handling

Store the product at -20°C or below upon receipt. The lyophilized probe is stable for >1 week at room temperature and for >12 months at -20°C. Reconstitute SPY555-Golgi using anhydrous DMSO. We recommend using newly or freshly opened and anhydrous DMSO to prepare the 1000x stock solution. In contact with air and moisture, DMSO produces decay products which can strongly reduce the shelf life of the probe in solution, even at -20°C. Keep the 1000x stock solution of the probe below -20°C after use. Vials should be allowed to warm to room temperature before opening. When reconstituted and stored properly, the 1000x stock solution is stable for 3 months. Note: DMSO solutions should be handled with particular caution as DMSO is known to facilitate the entry of organic molecules into tissues. Dispose of these reagents in compliance with all pertaining local regulations.

### Labelling Protocol

#### Important Note:

This protocol was optimized and validated using HeLa, MDCK, MCF-7, or RPE-1 cells adhering to coated glass or polymer dishes. For other cell lines, recommendations in this protocol should be used as a starting point, and optimal labeling conditions for each cell type should be determined empirically. High concentrations/long incubation times result in saturation of the golgi and leakage of the fluorescent signal into ER/endomembranes. Reducing the probe concentration or staining time by half usually restores a golgi specific signal.

**1. Prepare 1000x stock solution.** Add 50 µL of anhydrous DMSO to the SPY555-Golgi vial to prepare the 1000x stock solution (DMSO stock concentration is 300 µM). We recommend using newly or freshly opened and anhydrous DMSO to prepare the 1000x stock solution. In contact with air and moisture, DMSO produces decay products which can strongly reduce the shelf life of the probe in solution, even at -20°C. At this stage, the solution can be colored or not, this has no influence on the performance of the probe. After use, this solution should be stored at -20°C or below. Do not divide the 1000x stock solution into small aliquots, they will decay faster and the probe is not altered by multiple freeze-thaw cycles. When stored properly, this stock solution is stable for 3 months.

**2. Prepare the staining solution.** Dilute SPY555-Golgi to 1x (i.e. a final concentration 300 nM) in your usual cell culture medium (e.g. DMEM + 10% fetal bovine serum) and vortex briefly. If the dilution is not performed in a single step, please use DMSO to prepare the intermediate dilution as using aqueous buffers to prepare the intermediate dilution may lead to the formation of probe aggregates. Proceed quickly to step 3. Since staining efficiency can depend on the cell line, it is recommended to stain cells with 1000x dilution at the first attempt and then optimize the SPY555-Golgi dilution factor in further experiments until an optimal staining is achieved. Use only freshly made staining solution, and do not use it multiple times.

### **3. Cell preparation and staining.**

**3.1 Live cell staining:** Grow cells on coverslips, glass bottom dish or glass bottom multi-well plates as usual. When cells have reached the desired density, replace the culture medium by the **staining solution** freshly prepared under step 2 ensuring that all the cells are covered with the solution. Place the cells in the incubator at 37°C in a humidified atmosphere containing 5% CO<sub>2</sub> for 15-30 minutes\*\*. Perform a simple washing step consisting of replacing once the labelling solution by fresh culture medium which does not contain the probe.

**3.2 Fixed cell staining:** Grow cells on coverslips, glass bottom dish or glass bottom multi-well plates as usual. When cells have reached the desired density, the cells must be stained first with SPY555-Golgi as described in in step 3.1 and then fixed with 4% formaldehyde or 2% glutaraldehyde using standard protocols. **Post fixation staining cannot be performed**, SPY555-Golgi needs to be metabolically palmitoylated within the live cells to properly localize to the golgi.

**4. Cell imaging.** Imaging of SPY555-Golgi is best performed using standard Cy3/TMR settings. **Ca 3-4 hours after the staining, the signal intensity drops by 50%. Therefore, the probe should not be used for timelapse experiments exceeding 4h.**

\* Based on the following conditions: 0.5 ml staining solution / staining experiment with 1x probe concentration. The number of staining experiments can be further increased by reducing volume or probe concentration.

\*\* The recommended labelling time was determined for 2D cultured cells and may differ slightly depending on the cell line used or the sample type (e.g. spheroids, organoids or tissue).

Spirochrome products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply. Spirochrome products and product applications are covered by patents and patents pending. SPY is a registered trademark

**Limited Use Label License:** For research use only. Not intended for any animal or human therapeutic or diagnostic use. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Spirochrome will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Spirochrome is willing to accept return of the unused product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Spirochrome: Spirochrome AG, Chalberweidstrasse 4, 8620 Stein am Rhein, Switzerland, Email: info@spirochrome.com