

Product information: BD566-CA (SC220)


BD566-chloroalkane ligand for HaloTag™ labeling

Introduction


Product Description

BD566-CA (also named BD566-HTL) is a ChloroAlkane (CA) derivative of the bridged-bicycle-strengthened orange fluorophore BD566. The unique BriDye™ bicyclic ether bridge technology confers to BD566-CA an exceptional photostability. Its absorption and emission maxima are 566 nm and 594 nm, respectively with an extinction coefficient of 108 000 M⁻¹ cm⁻¹ and a quantum yield of 0.85 (HT7 bound values). BD566 is highly cell-permeable, highly photostable, and compatible with standard Cy3 or TMR filter settings for widefield, confocal, SIM, or STED imaging in living or fixed specimens.

ChloroAlkane (CA) is the substrate of the self-labeling tag HaloTag™. Upon reaction with a CA derivative, HaloTag™ forms a covalent bond with the substrate, enabling permanent attachment of a fluorescent label to any protein of interest (POI) expressed as a HaloTag™ fusion. Contains 1 vial of BD566-CA (10 nmol, lyophilized).



BD566-CA structure



BD566-CA absorption (solid) and emission (dashed) spectra when bound to HaloTag 7

Properties

Absorbance maximum λ_{abs}	566 nm
Fluorescence maximum λ_{fl}	594 nm
Works on fixed cells?	yes
Quantity	10 nmol
Fluorescence lifetime (free dye)	3.2 ns
MW	772.3 g/mol
Shipping	room temperature
Storage	-20°C

Storage & Handling

Store the BD566-CA at -20°C or below upon receipt. The lyophilized HaloTag™ substrate is stable for >1 week at room temperature and for at least 6 months at -20°C. Reconstitute BD566-CA using anhydrous DMSO. We recommend using newly or freshly opened and anhydrous DMSO to prepare the stock solution. In contact with air and moisture, DMSO produces decay products which can strongly reduce the shelf life of the compound in solution, even at -20°C. Keep the stock solution of the BD566-CA below -20°C after use. Vials should be allowed to warm to room temperature before opening. When reconstituted and stored properly, the stock solution is stable for at least 3 months. Note: DMSO solutions should be handled with particular caution as DMSO is known to facilitate the entry of organic molecules into tissues. Dispose of these reagents in compliance with all pertaining local regulations.

Labelling Protocol

Note: Recommendations in this protocol should be used as a starting point, and optimal labeling conditions for each cell type should be determined empirically.

1. Prepare DMSO stock solution. Add 50 µL of anhydrous DMSO to the BD566-CA vial to prepare a 1000x (200 µM) stock solution. We recommend using newly or freshly opened and anhydrous DMSO to prepare the DMSO stock solution. In contact with air and moisture, DMSO produces decay products which can strongly reduce the shelf life of the substrate in solution, even at -20°C. At this stage, the solution can be colored or not, this has no influence on the performance of BD566-CA. After use, this solution should be stored at -20°C or below. Do not divide the DMSO stock solution into small aliquots, they will decay faster and the compound is not altered by multiple freeze-thaw cycles. When stored properly, this stock solution is stable for 3 months.

2. Prepare the staining solution. Dilute the BD566-CA DMSO stock solution 1:1000 (final concentration 200 nM) in your usual cell culture medium (e.g. DMEM + 10% fetal bovine serum) and vortex briefly. Proceed quickly to step 3. If the dilution is not performed in a single step, use DMSO to prepare the intermediate dilution as using aqueous buffers to prepare the intermediate dilution may lead to the formation of aggregates. Since staining efficiency can depend on the cell line, it is recommended to stain cells with 1:1000 dilution at the first attempt and then optimize the dilution factor in further experiments until an optimal staining is achieved. The usual concentration range for live cell labelling is 100-500 nM. Use only freshly made staining solution and do not use it multiple times.

3. Cell preparation and staining. Grow cells transiently transfected or stably expressing a Halotag™* fusion-protein on coverslips, glass bottom dish or glass bottom multi-well plates as usual. When cells have reached the desired density or expression level of the Halotag™* fusion-protein, replace the culture medium by the **staining solution** freshly prepared under step 2 ensuring that all the cells are covered with the solution. Place the cells in the incubator at 37°C in a humidified atmosphere containing 5% CO₂ for min 30 minutes, although 1h is recommended.

Note: Before imaging, BD566-CA stained cells can be fixed by any fixation method after the labelling step is completed. Additional immunolabeling or probe labeling can be performed after the fixation step using standard protocols.

4. Cell imaging. Imaging of BD566-CA labeled cells is best performed using standard Cy3 or TMR settings. After labelling the live cells, a simple washing step consisting of replacing once or twice the labelling solution by fresh culture medium which does not contain the substrate will suppress the background signal from unreacted ligand. The staining will last depending on your Halotag™* fusion protein stability and turnover rate.

*Halotag™ is a registered trademark of Promega

BD566-CA and BD626-CA are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply. BD566-CA and BD626-CA products and product applications are covered by patents and patents pending. BD566-CA and BD626-CA products are proposed in partnership with Nanjing Genvivo Biotech Co. Ltd. products under an exclusive license. BriDye™ is a registered trademark of Nanjing Genvivo Biotech Co. Ltd.

Limited Use Label License: For research use only. Not intended for any animal or human therapeutic or diagnostic use. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Spirochrome will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Spirochrome is willing to accept return of the unused product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Spirochrome: Spirochrome AG, Chalberwidstrasse 4, 8620 Stein am Rhein, Switzerland, Email: info@spirochrome.com